



# LUXATING HYBRIDS

PRECISION AND DURABILITY COMBINED

## FACT

The use of unsuitable surgical instruments for tooth extractions can lead to undesirable results and patient dissatisfaction.

## CHALLENGE

Extracting teeth or broken root tips without causing trauma to the buccal plate can become difficult without the correct instruments.

## SOLUTION

Luxating Hybrids were designed with thin, precisely angled tips to aid clinicians in preserving the buccal plate during atraumatic extractions.



INCREASED INSTRUMENT CONTROL WITH ERGONOMIC, SATIN STEEL® HANDLE



## POINTS OF PERFORMANCE

- **Improved apical and interproximal access** with thin working tips
- **Easily severs the periodontal ligament** with finely honed, sharp cutting edges
- **Enhanced oral adaptation** with 7 precisely manufactured tip options
- **Longevity and corrosion resistance** afforded by hand-crafting from the finest grade surgical stainless steel

SEE FULL LISTING OF LUXATING HYBRID INSTRUMENTS ON THE BACK.

How the best perform

**Hu-Friedy**



## LUXATING HYBRIDS

**STRAIGHT SPADE (EPTSSP)**



**SMALL STRAIGHT (EPTSMS)**



**SMALL CURVED (EPTSMC)**



**SMALL DISTAL (EPTSMD)**



**SMALL MESIAL (EPTSMM)**



**LARGE DISTAL (EPTLGD)**



**LARGE MESIAL (EPTLGM)**



For clinicians who are interested in the entire Luxating Hybrid line, Hu-Friedy has created a kit with all seven instruments for atraumatic extractions.

It also includes a sharpening stone to keep these instruments performing like new.



### LUXATING HYBRID KIT

### EPTKIT

Small Signature Series® 8 Instrument Cassette, Red	IM5080
Cylindrical Ceramic Fine Grit Sharpening Stone	SSKC
Luxating Hybrid – Straight Spade	EPTSSP
Luxating Hybrid – Small Straight	EPTSMS
Luxating Hybrid – Small Curved	EPTSMC
Luxating Hybrid – Small Distal	EPTSMD
Luxating Hybrid – Small Mesial	EPTSMM
Luxating Hybrid – Large Distal	EPTLGD
Luxating Hybrid – Large Mesial	EPTLGM

Visit us online at [Hu-Friedy.com](http://Hu-Friedy.com)

©2012 Hu-Friedy Mfg. Co., LLC. All rights reserved. M4656/0812

**Hu-Friedy**